



## FREEZING COLOSTRUM

### ▶▶ CRITERIA:

- 1<sup>st</sup> milking
- "Excellent" colostrum
- Freshly collected

### ▶▶ STEPS:

- Large Ziploc freezer bags
- Labelled (cow, date, quality)
- Approx. 1 L per bag
- Freeze flat
- Store for maximum of 12 months

## THAWING COLOSTRUM

### ▶▶ STEPS:

- Put the Ziploc in a second bag (to prevent leakage)
- Place Ziploc in a bucket of tepid water (Max 40 °C) (water bath)
- Let thaw (15-20 min.)
- Transfer into a feeding bottle
- Heat to 38 °C
- Feed to the calf